Acupressuurmat en kussen

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
ACUPRESSURE MAT THERAPY

The acupressure mat with pillow is equipped with over 8,000 pressure points (spikes) that can stimulate the skin and muscles along meridian lines throughout your body. This has many pleasant advantages, including the:

- Promote general vitality
- Relaxing body and mind
- Stimulate the production of endorphins
- Promote blood circulation
- Reducing pain and muscle tension
- Reducing stress
- Increasing energy
- Improving sleep


For best results, we recommend lying on the mat with bare skin. If you don't like this, you can also wear thin clothing or put a thin fabric on the mat to get your body used to the feeling of acupressure. Be sure to lie on and get off your mat carefully, as the pressure points are very sharp.

Give in to the acupressure experience

You may find the feel of the pressure points uncomfortable at first because you still need to get used to the sensation of the points. Take a deep breath and relax and within a few minutes, this uncomfortable feeling will disappear on its own, after which you will be able to truly experience the benefits of acupressure.

Tips for use

These are just some of the basic exercises you can use the mat for, but the possibilities are endless. Experiment and research yourself. Focus mainly on an area you want to work with, surrender and experience! When you first use the mat on a new area of your body, this area will often be sensitive. The body is intelligent; It will get used to the feeling over time and surrender to it.


Research what works for your body. For some areas of the body, just a few minutes on the mat is enough, while other areas, such as the back, may require more time. The areas of the body that hurt the most while using the mat usually need this use the most. The pain may reflect tension stored in these areas.

The basic exercises

The acupressure mat and pillow can be used for many exercises and in many postures. Here are some basic exercises suitable for both those with and those without experience.


Feet

- Effect: relieves tired feet, makes you feel more energetic and brings a sense of recovery. For a great start and end to your day!
- Instruction: Stand on the mat or sit on a chair with your feet on the pillow. Wear (thin) socks in the beginning and slowly build up to bare feet.
- Time: 5-10 minutes


Neck and back

- Effect: can reduce back pain neck and shoulder tension and migraines. Reduces stress, promotes relaxation and sleep.
- Instructions: Lie on a flat surface with your back on the mat and the pillow directly under your neck. For higher intensity in the lower back, you can bend your knees and put your feet flat on the ground. Fine as part of your evening ritual!
- Time: Beginners: 10-20 minutes. Advanced: up to 45 minutes


Chest and abdomen

- Effect: Brings recovery after workouts, supports digestion and promotes circulation in respiratory, abdominal and intestinal muscles.
- Instructions: Lie on your stomach and rest your head off the mat on your hands. You can also lie with one of your cheeks on the mat or on the accompanying pillow to relax your face and jaw at the same time.
- Time: 10-20 minutes


Lower back and gluteal muscles

- Effect: Relieves tension in the lower back, provides recovery after lower back and leg workouts. Promotes good posture, useful for office use!
- Instructions: Place the mat on a chair with a slight fold so that both gluteus and lower back make contact with the mat. Sit gently against the mat and your chairrest.
- Time: 10-20 minutes

Legs and hamstrings

- Effect: Targets various pains or tensions in the lower body and works to promote recovery after workouts or an intense work day. Supports healthy circulation, calms the muscles and provides targeted relaxation.
- Instructions: Sit on the edge of the mat and stretch your legs out in front of you. To work on your flexibility and stimulate circulation, reach for your shins or toes and maintain this for 20-30 seconds per minute each time.
- Time: 10-15 minutes

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